

FAQs

WARNING: If the instruments with the branches are placed in the modules (on the integrated magnets) for a long time, the instruments may become magnetized. Please make sure that you either remove the instruments after each exercise or place them in the box so that there is a sufficient distance from the modules. LTB Germany Ltd. assume no liability for any magnetization.

■ How much space is needed for the Luebeck Toolbox setup?

The required area for a practicing person and a Toolbox together with the monitor is about 2x1m. The Toolbox can be positioned on any flat surface with the dimension 40x50cm (you may have to add some kind of elevation).

■ Is the Luebeck Toolbox transportable?

The Toolbox can be comfortably carried by one Person and set up in different locations.

■ How long does it take to Setup the Luebeck Toolbox?

If all the required components are available, the installation of a workplace takes no longer than 15 minutes.

■ Is the monitor included?

The monitor (conventional TV monitor with HDMI connection) is not included.

■ What do I need in addition to the Toolbox for practicing?

One monitor, two trocars, the instruments (one curved scissors, two Overholt graspers, two needle drivers) and the sutures (braided, 3-0, SH needle).

■ Can I get certified upon successful completion of the curriculum?

The curriculum is accredited with 24 CME points. Find out more at http://www.luebeck-toolbox.com/files/pdf/LTB_Fortbildungspunkte_DE.pdf

■ What options do I have to contact LTB Germany Ltd.?

You can contact us by email at the following addresses:

ltb@luebeck-toolbox.com or support@luebeck-toolbox.com

■ Where can I find the results of the studies that lead to the conception of the Luebeck Toolbox curriculum?

At Pubmed you will find all previously published original works on the Lübeck Toolbox.

■ Can I do other exercises with the Luebeck Toolbox than those included in the curriculum?

In principle, yes. The modules are fixed with four small magnets at corresponding points of the bottom plate of the Toolbox. If the user manages to fixate other training modules on the base plate, they of course could also be used. However, the training success of the Luebeck Toolbox Curriculum is solely validated for the given exercises. If a suitable plastic container is used, even organ procedures could be performed. However, the LTB Germany Ltd does not carry liability for the required ethical and hygienical standards.

■ Are the Toolbox and the modules readily assembled?

All you have to do is install the camera in the Toolbox and place the modules for each exercise with the delivered materials (plastic pegs, rubber bands, non-woven fabric, etc). Before insertion of the trocars, the rubber membranes must be cut.

■ How long does the shipping take?

Maximum 4 weeks, on average we expect 3 weeks. If required, the delivery can also be faster using express shipping. A tracking number will be provided.

■ How long is the warranty?

The warranty for the Luebeck Toolbox including the camera and the 4 modules is one year. Further details can be found in under terms and conditions.

■ Can I withdraw from the purchase?

The buyer may withdraw from the purchase within two weeks without stating any reasons. Further details can be found under terms and conditions.

■ Where is the Luebeck Toolbox manufactured?

The Toolbox is manufactured and assembled by a renowned company in southern Germany. The modules are manufactured near Hong Kong by established manufacturers. A thorough quality control of all components is carried out in a functional test of the entire system before shipping.